

## RNA Extraction Solution

**Product #:** G3013-100ML

**Quantity:** 100 mL

### Product Description

- RNA extraction solution is a ready-to-use reagent, designed to isolate high quality total RNA from cell and tissue samples.
- This product utilizes a monophasic solution of phenol, guanidine isothiocyanate, and other proprietary components which facilitate the isolation of a variety of RNA species of all sizes (miRNA to large RNA products).
- RNA extraction solution maintains RNA integrity due to highly effective inhibition of RNase activity while disrupting cells and dissolving cell components during sample homogenization.
- Purified RNA products can be directly used for subsequent molecular biology applications.

### Storage

- Store at 4°C in a dark environment for up to 12 months.

### Components Provided

RNA Extraction Solution      100 mL

### User Supplied Components

---

Nuclease Free H<sub>2</sub>O

---

Chloroform

---

Isopropanol

---

75% Ethanol

---

Centrifuge and rotor capable of reaching 12,000 x *g* and 4 °C

---

**Important Preparatory Notes (PLEASE READ CAREFULLY)**

1. Perform all steps at room temperature (20–25°C) unless otherwise noted.
2. Use cold **RNA Extraction Solution** if the starting material contains high levels of RNase, such as spleen or pancreas samples.
3. Use disposable, individually wrapped, sterile plasticware and sterile, disposable RNase-free pipettes, pipette tips, and tubes.
4. Wear disposable gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin; change gloves frequently, particularly as the protocol progresses.
5. Always use proper microbiological aseptic techniques when working with RNA.

**Samplly Lysis****Tissues:**

1. Add 1 mL of **RNA Extraction Solution** per **50 – 100 mg** of tissue.
2. Homogenize sample until there are no visible tissue clumps.
3. Centrifuge at **12,000 x g** for **10 min** at **4°C** to separate the tissue fragments.
4. Transfer the supernatant to a new centrifuge tube.

**Adherent Cells:**

1. Remove growth media and wash cells with **ice-cold 1X PBS**.
2. Remove PBS wash thoroughly.
3. Add 0.25 to 1 mL of **RNA Extraction Solution** per **1 × 10<sup>5</sup>–10<sup>7</sup> cells** directly to the culture dish to lyse the cells.
4. Pipet the lysate up and down several times to homogenize.
5. Transfer the sample to a new centrifuge tube.
6. Let the sample stand at **room temperature** for **5 min** before proceeding.

**Cell Suspension:**

1. Pellet the cells by centrifugation and discard the supernatant thoroughly.
2. Add **0.25 to 1 mL** of **RNA Extraction Solution** per **1 × 10<sup>5</sup>–10<sup>7</sup> cells**.
3. Pipet the lysate up and down several times to homogenize.
4. Transfer the sample to a new centrifuge tube.
5. Let the sample stand at **room temperature** for **5 min** before proceeding.

**RNA Purification Protocol**

1. Add **380  $\mu$ L** of **chloroform** per **1 mL** of **RNA Extraction Solution**.
2. Vortex vigorously for 15 s, and let sample stand at **room temperature** for **3 min**.
3. Centrifuge sample at **12,000 x g** for **15 min** at **4°C**.
4. Carefully transfer the upper colorless aqueous phase to a new centrifuge tube.
  - **Note:** *Approximately 500-550 $\mu$ L of aqueous solution can be extracted per 1 mL of RNA Extraction Solution.*
  - **Note:** *the interphase and organic phases can be saved for DNA or protein isolation.*
5. Precipitate RNA by adding **550 $\mu$ L** of **isopropanol**. Gently invert several times.
6. Place tube at **-20°C** for **15 minutes** to enhance the precipitation of the RNA.
7. Centrifuge at **12,000 x g** for **10 min** at **4°C**.
8. Total RNA precipitate forms a white gel-like pellet at the bottom of the tube. Discard the supernatant using a micropipettor.
9. Wash the pellet by adding **1 mL** of **75% ethanol** and mix by inversion.
10. Centrifuge at **12,000 x g** for **5 min** at **4°C**.
11. Discard the supernatant using a micropipettor.
12. Briefly centrifuge (**5,000 x g** for **3-5 s**), and carefully remove the remaining liquid using a micropipettor.
13. **Air-dry** the RNA pellet by opening the centrifuge tube for **5 min**.
14. After the RNA pellet has slightly dried, add 20 – 50  $\mu$ L of nuclease-free water to fully dissolve the RNA. Mix by pipetting or vortex.
  - **Note:** *Be careful not to let the RNA dry too much, or the pellet is difficult to dissolve.*
  - **Optional:** *Incubate on a thermomixer for 10 min at 1000 rpm and 37°C.*
15. Determine the RNA yield and purity.
16. Store the purified RNA at **-80°C** for long-term storage. Avoid repeated freeze-thaw of the RNA.

**NOTES:**

1. **DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
2. **Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
3. ©2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
4. For Research Use Only. Not for use in diagnostic procedures.