## **USER GUIDE** G3640-50T

### **RNA Extraction Kit (Spin Column Based)**

Product #: G3640-50T Quantity: 50 reactions

#### **Product Description**

- Designed for rapid and cost-effectives purification of RNA from cultured cells (adherent and non-adherent) and animal tissues.
- High purity RNA extraction is achieved by combining proprietary silica-based membrane technology in the form of a convenient spin column. The kit utilizes two columns: gDNA Eraser Spin Column (to remove genomic DNA) and an RNA Spin Column (to bind RNA).
- The procedure avoids phenol and chloroform. The purified RNA contains no protein and no genomic DNA contamination.
- The high-quality RNA can be used in Northern blotting, dot blotting, mRNA purification, *in vitro* translation, RNase protection assays, RT-PCR, Real Time RT-PCR, cDNA library construction and other molecular experiments including direct sequencing.

Buffer RL1	30 mL	
Buffer RW1	12 mL (note: add 18 mL 100% ethanol per bottle before use)	
Buffer RW2	20 mL (note: add 80 mL 100% ethanol per bottle before use)	
Nuclease-free Water	12 mL	
gDNA Eraser Spin Column	50 columns	
RNA Spin Column	50 columns	
10×DNase Buffer	500 μL	

#### Kit Components – BAG #1 (store at room temperature)

#### Kit Components – BAG #2 (store at -20°C)

DNase	250 μL
DTT Solution	1.2 mL

100% Ethanol	
Isopropanol	

## **USER GUIDE** G3640-50T

#### **Important Preparatory Notes (PLEASE READ CAREFULLY)**

- 1. If precipitate forms in **Buffer RL1** during storage, solution should be redissolved by incubating at 37°C, then cooled back down to 25°C before use.
- 2. Add **18 mL 100% ethanol** to the **Buffer RW1** and mark the bottle on the label.
- 3. Add **80 mL 100% ethanol** to the **Buffer RW2** and mark the bottle on the label.
- 4. Add appropriate volume of **DTT Solution** into **Buffer RL1** to a final concentration of 4% (i.e. 40  $\mu$ L DTT Solution per 1 mL of Buffer RL1). This addition should be prepared immediately before use. This mixture can be stored at 4°C for one month.
- 5. Always keep buffer bottles tightly closed and keep columns sealed in the enclosed zip-lock bag.
- 6. Proper laboratory safety practices should be employed, including the use of lab coats, gloves, and eye protection.

#### **RNA Extraction Protocol**

- 1. Sample preparation:
  - Non-adherent suspension cells: Transfer cells to an appropriately sized RNase-free tube and centrifuge at 1,000 rpm for 5 minutes at 4°C to pellet. Discard the growth medium from the tube. Add 500 μL of Buffer RL1 prepared with DTT solution to the sample. Pipet up and down until no significant pellet. Let the sample stand at room temperature for 2 minutes.
  - Adherent cells: Remove the growth medium from the cells and wash cells with 1×PBS buffer (pH 7.4). Add 500 µL of Buffer RL1 prepared with DTT solution evenly over the cell monolayer. Pipet the cells up and down several times. Transfer the lysate to a 1.5 ml RNase-free tube. Let the sample stand at room temperature for 2 minutes.
  - Animal tissue: Transfer 2-20 mg of fresh or cryopreserved animal tissue to an appropriately sized RNase-free tube. Add two or three 4 mm stainless

# **WENDER GUIDE CORP. G3640-50T**

steel grinding beads (Cat# G0104) per sample. Immediately add **500**  $\mu$ L of **Buffer RL1** prepared with **DTT solution** to the sample. Use an automated homogenizer (TissueLyser) to homogenize thoroughly (2 to 5 min @ 20-50 Hz). Let the sample stand at room temperature for 2 minutes.

- 2. Transfer sample to **gDNA Eraser Spin Column**. Centrifuge at 12,000 rpm at 4°C for 1 minute.
- 3. Discard the **gDNA Eraser Spin Column** while keeping the flow-through in the collection tube.
- 4. Add **2/3 volume** of **isopropanol** (*note: visible precipitate may form after adding isopropanol*) to the flow-through in the collection tube. Pipet up and down to mix well.
- 5. Transfer up to 600  $\mu\text{L}$  of the sample (including any remaining precipitate) to the RNA Spin Column.
- 6. Centrifuge at 12,000 rpm for 30 seconds at room temperature. Discard the flow-through, and reinsert the **RNA Spin Column** in the same collection tube.
- Add 500 μL Buffer RW1 to the spin cartridge. Centrifuge at 12,000 rpm for 30 seconds at room temperature, and discard the flow-through. Reinsert the RNA Spin Column in the same collection tube.
- Add 600 μL Buffer RW2 to the spin cartridge (*note: add Buffer RW2 along the wall of the RNA Spin Column to help flush the residual salt*). Centrifuge at 12,000 rpm for 30 seconds at room temperature, and discard the flow-through. Reinsert the RNA Spin Column in the same collection tube.
- 9. Optional DNase digestion: The **gDNA Eraser Spin Column** can effectively remove most of the genomic DNA. For some tissues with high DNA content (such as liver, kidney, spleen), perform the optional DNase digestion steps:
  - Add **5**  $\mu$ L DNase to **5**  $\mu$ L **10X** DNase Buffer and **40**  $\mu$ L RNase-free Water for a total of 50  $\mu$ L in a new 1.5 mL RNase-free centrifuge tube. Mix by gently pipetting up and down several times.
  - Pipet all 50 μL of the DNase mixture onto the center of the RNA Spin Column membrane. Incubate at room temperature for 15 minutes.



- Add 500 μL Buffer RW2 to the spin cartridge. Centrifuge at 12,000 rpm for 30 seconds at room temperature, and discard the flow-through. Reinsert the RNA Spin Column in the same collection tube.
- 11. Repeat Step 10.
- 12. Centrifuge the **RNA Spin Column** at 12,000 rpm for 2 minutes to dry the membrane with RNA.
- 13. Discard the flow-through and the Collection Tube and insert the **RNA Spin Column** into a **new 1.5 mL RNase-free centrifuge tube**. Incubate at room temperature for 3-5 minutes to completely evaporate residual ethanol in the RNA Spin Column.
- 14. Add **50-100 μL Nuclease-Free Water** to the center of the **RNA Spin Column**. Incubate at room temperature for 5 minutes.
- 15. Centrifuge the **RNA Spin Column** for 2 minutes at 12,000 rpm at room temperature to elute the RNA from the membrane into the RNase-free centrifuge tube.
  - If more yield is required, the elution can be re-applied into the center of the membrane. Incubate for 5 minutes at room temperature and centrifuge at 12,000 rpm for 2 minutes to elute the RNA.
- 16. Store the purified RNA at -80°C for long-term storage. Avoid repeated freezing and thawing of the DNA.

## **WER GUIDE** G3640-50T

#### **Approximate RNA Yield:**

Sample species	Sample names	RNA yield
Animal tissues	Mouse liver	30-45 µg/10 mg
	Mouse spleen	20-30 µg/10 mg
	Mouse lung	10-20 µg/10 mg
	Mouse kidney	20-30 µg/10 mg
	Mouse heart	5-10 µg/10 mg
	Mouse thymus	10-20 µg/10 mg
	Mouse brain	5-10 μg/10 mg
	Mouse pancreas	5-15 μg/10 mg
Cells	Hela	8-15 µg/10 <sup>6</sup> cells

NOTES:

- 1. DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
- 2. Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
- 3. © 2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
- 4. For Research Use Only. Not for use in diagnostic procedures.