

Recombinant Anti β -Actin Mouse Monoclonal Antibody

Product #: Bi2M-GB15001

Quantity: 1 x 25 μ L

Product Overview

Target	beta-Actin (β -Actin)
Immunogen	KLH conjugated synthetic peptide corresponding to mouse β -Actin
Species Reactivity	Human, Mouse, Rat
Host/Isotype	IgG1
Class	Monoclonal
Type	Primary Antibody
Conjugate	None
Form	Liquid
Concentration	0.01 mg/mL
Purification	Antigen affinity chromatography
Storage Buffer	PBS with 0.02% sodium azide, 100 μ g/mL BSA and 50% glycerol
Contains	No preservative
Applications	Western Blot, Immunohistochemistry

Product Description

Affinity purified mouse monoclonal (heavy and light chain) antibody targeting β -Actin.

Storage

Store at -20°C for long-term storage. Once thawed, store at 4°C . Do not freeze-thaw.

Product Usage Information - Recommended Dilutions

Western Blot	1:500 – 1:1,000
Immunohistochemistry	1:50 – 1:200

NOTES:

- DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
- Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
- ©2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
- For Research Use Only. Not for use in diagnostic procedures.