

Goat Anti-Rabbit Secondary Antibody – HRP Conjugated

Product #: Bi2M-GB23303

Quantity: 1 x 100 µL

Product Overview

Species Reactivity	Rabbit
Host/Isotype	Goat/IgG (H+L)
Class	Polyclonal
Type	Secondary Antibody
Conjugate	HRP – Horseradish Peroxidase
Form	Liquid
Concentration	0.01 mg/mL
Purification	Antigen affinity chromatography
Storage Buffer	PBS, pH 7.6
Contains	No preservative
Applications	Western Blot, Immunohistochemistry

Product Description

Affinity purified goat anti-rabbit IgG (heavy and light chain) antibody conjugated to horseradish peroxidase (HRP) for chemiluminescent detection. This product has been verified for use in Western blot and immunohistochemistry-paraffin applications. This antibody may cross-react with immunoglobulins from other species.

Storage

Store at -20°C for long-term storage. Once thawed, store at 4°C. Do not freeze-thaw.

Product Usage Information - Recommended Dilutions

Western Blot	1:2,000 – 1:10,000
Immunohistochemistry	1:200 – 2,000

NOTES:

- DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
- Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
- ©2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
- For Research Use Only. Not for use in diagnostic procedures.