

dNTP mix (10mM)

Product #: Bi2M-dNTP

Quantity: 1 mL

Product Description

- This product is a mixture of high-purity dATP, dTTP, dGTP, and dCTP, each at a concentration of 10 mM.
- Ready for use in PCR, RT-PCR, real-time PCR, DNA sequencing, DNA labeling, and other molecular biology experiments.
- Packaged from imported high-purity reagents with HPLC-detected purity over 99%, free of DNase and RNase.
- For a standard 50 μ L PCR reaction system, the recommended dosage is 1 μ L (final concentration of 200 μ M).

Storage

- Store at -20°C.

NOTES:

1. **DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
2. **Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
3. ©2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
4. For Research Use Only. Not for use in diagnostic procedures.