



USER GUIDE

G4700-500ML

Molecular Water (DNase/RNase/Protease Free)

Product #: G4700-500ML

Quantity: 500 mL

Product Description

- Nuclease and protease free ultra pure water.
- Non-DEPC treated and deionized.
- Manufactured using multiple distillation followed by 0.1 μ m filtration and sterilization.
- Can be used for downstream nucleic acid related molecular biology experiments and routine biology experiments.

Storage

- Room temperature

NOTES:

1. **DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
2. **Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
3. ©2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
4. For Research Use Only. Not for use in diagnostic procedures.