

Bovine Serum Albumin (BSA)

Product #: G5001-100G

Quantity: 100g

Product Description

- Bovine Serum Albumin (BSA) is a protein derived from cow's blood serum.
- It is commonly used in various biochemical applications due to its high purity and stability.
- BSA serves as a standard protein in assays, acting as a carrier protein for hormones, enzymes, and fatty acids.
- It is utilized in immunoassays, protein quantification, enzyme assays, and cell culture applications.
- BSA helps stabilize enzymes and antibodies, preventing their denaturation and maintaining their activity.

Storage

- Store at 2-8°C and protect from light.

NOTES:

1. **DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, MEDIRES CORP. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.
2. **Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable limited Use Label Licenses.
3. ©2024 MEDIRES CORP. All rights reserved. All trademarks are the property of MEDIRES CORP. and its subsidiaries.
4. For Research Use Only. Not for use in diagnostic procedures.